



How to clear/clean Starkey Halo Hearing Aid microphone screens.

How to clear/clean Starkey Halo Hearing Aid microphone screens.

Written By: Jon Hirsch



INTRODUCTION

After a period of wearing the hearing aids, the microphone stops working and all you can hear is the bluetooth input. This video shows how to remove and clean the cover over the mics.

Remove battery

Remove in-ear speaker wire

Remove cover

Clean cover

Replace all parts

[video: <https://youtu.be/fWyNUHK2Wg8>]

To reassemble your device, follow these instructions in reverse order.